

Native Human-Derived ECM (Lyophilized)

Cat. No. MT005

Lot. No. (See product label)

Product Description	Due to its human origin and minimal processing, Native Human- Derived ECM retains the native matrix proteins and growth factors found in human tissue including collagen, elastin, laminin, glycosaminoglycans and many other matrix proteins.
Application	Native Human-Derived ECM is ideal for many applications including coating tissue culture surfaces to support cell attachment and growth or for making hydrogels.
State	Lyophilized
Biological Source	Human
Storage/Stability	Stored at -20 °C upon receipt
Shelf Life	6 months after opening
Note	Making 3D Hydrogels Note: All procedures should be performed in a cold room or on ice. All reagents should be cooled prior to making the hydrogel. To prevent premature gelation, keep solutions between 2-10°C. Concentrations will need to be tested to determine the optimal gel concentration for each application. Native Human-Derived ECM has been tested and found to form a gel at ≥ 3 mg/mL final concentration.
	1. Reconstitute Native Human-Derived ECM(Lyophilized) in 10-20 mM HCl with a pH of 1.9-2.1. It is recommended to use ultrapure sterile water to make the 10-20 mM HCl.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

www.matexcel.com/hydrogel

Address: SUITE 210, 17 Ramsey Road, Shirley, NY 11967

Email: info@matexcel.com Tel: 1-631-869-4956 Fax: 1-631-910-2166



- 2. Prepare 10x culture medium or 10x PBS with antibiotic and antimycotic.
- 3. Combine 9 volumes of chilled solution with 1 volume of chilled 10x concentrated medium or 10x PBS.
- 4. Gently swirl or pipette the mixture repeatedly.
- 5. Adjust the pH to 7.0-8.0, preferably around 7.4. Concentration and pH will affect the speed and strength of the final gel. Once pH has been adjusted, cells can be added to the solution if desired.
- 6. Incubate at 37°C for 30-120 minutes for gel formation depending on final concentration.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

www.matexcel.com/hydrogel

Address: SUITE 210, 17 Ramsey Road, Shirley, NY 11967

Email: info@matexcel.com Tel: 1-631-869-4956 Fax: 1-631-910-2166