



SG-Dextran

Cat. No. PS012

Lot. No. (See product label)

Product Description

SG-Dextran contains dextran functionalized with thiol-reactive groups that can be used together with PEG-Link or CD-Link to set up slow gelling, cell-compatible hydrogels.

10x CB (pH 7.2) buffer allows control of pH and osmotic conditions during gel formation in the presence of cells.

Application

3-D cell cultures

Embedment and immobilization of cells or tissues

Subsequent recovery of specimens by hydrogel dissolution

Components

SG-Dextran: 3 x 170 μ L, 30 mmol/L

10x CB (pH 7.2): 600 μ L

Water: 4 x 1500 μ L

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

www.matexcel.com/hydrogel

Address: SUITE 210, 17 Ramsey Road, Shirley, NY 11967

Email: info@matexcel.com Tel: 1-631-869-4956 Fax: 1-631-910-2166