



## CD-Link

Cat. No. MT010

Lot. No. (See product label)

<b>Product Description</b>	The CD-Link molecule consists of a polyethylene glycol-peptide conjugate carrying a thiol at each end.
<b>Application</b>	Thiol groups on CD-Link form stable thioether bonds with the thiol-reactive groups on the polymers, which results in the formation of the hydrogel.
<b>Components</b>	CD-Link, lyophilized: 200 $\mu\text{L}^*$ , 20 mmol/L $^*$ Water: 600 $\mu\text{L}$

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

[www.matexcel.com/hydrogel](http://www.matexcel.com/hydrogel)

Address: SUITE 210, 17 Ramsey Road, Shirley, NY 11967

Email: [info@matexcel.com](mailto:info@matexcel.com) Tel: 1-631-869-4956 Fax: 1-631-910-2166